

CASE MEDICAL, INC.

STERILIZATION CONTAINER SYSTEM VIDEO POST-TEST

1. Rigid containers should be completely disassembled and washed:
 - (a) When containers are used to transport contaminated instruments
 - (b) When containers are transported outside of the facility
 - (c) After each use
 - (d) Once a week

2. Cleaning anodized aluminum containers require the use of detergents that are:
 - (a) Acid
 - (b) Neutral pH
 - (c) Alkaline
 - (d) Highly Alkaline

3. Inspection should be performed to verify that gaskets in the lid and filter retention plate(s) are not cracked and are properly and firmly affixed:
 - (a) After each use
 - (b) Once a week
 - (c) On a monthly schedule
 - (d) During preventive maintenance

4. Internal chemical indicators should be placed inside each container prior to sterilization:
 - (a) In the center of the basket
 - (b) In the corner of the inside basket
 - (c) On left side of the basket only
 - (d) On right side of the basket only

5. When using sealed containers, additional exposure time may be required in this method of sterilization:
 - (a) Pre-vacuum steam
 - (b) EO
 - (c) Gas plasma
 - (d) Gravity displacement steam

CASE MEDICAL, INC.

STERILIZATION CONTAINER SYSTEM VIDEO POST-TEST

6. Rigid containers should be placed on the sterilizer cart:
 - (a) On the bottom shelf only
 - (b) On the top shelf only
 - (c) On their side
 - (d) In a flat position

7. Placement of tamper-evident seals/security locks on containers should be done:
 - (a) After the container is sterilized
 - (b) Prior to sterilization
 - (c) Just before distributing the container
 - (d) Whenever tape is not used to seal the container

8. Following the sterilization process, the cart should be removed from the sterilizer and placed in cool down in order to:
 - (a) Allow for placement of tamper evident seals/security locks on all containers
 - (b) Maximize the potential of condensate formation within the container
 - (c) Minimize the potential of condensate formation within the container
 - (d) Neutralize condensate formation within the container

9. At point of use before opening the container, always verify that the correct set has been selected and:
 - (a) Tamper-evident seals/locks are intact
 - (b) Filter retention plates are intact
 - (c) Gaskets are undamaged
 - (d) Acceptability of internal chemical indicator

10. Users of rigid container systems should become familiar with information contained in the:
 - (a) OSHA Standard
 - (b) JCAHO Standard
 - (c) EPA Standard
 - (d) AAMI Standards

CASE MEDICAL, INC.

STERILIZATION CONTAINER SYSTEM POST-TEST ANSWER SHEET

1. C
2. B
3. A
4. B
5. D
6. D
7. B
8. C
9. A
10. D