

Case Medical®

Square Polypropylene Sterilization Filters

The Case Medical polypropylene, non-cellulose disposable sterilization filter allows for penetration of the sterilizing agents (steam or gas) into the rigid container. It is hydrophobic, thereby reducing the risk of residual humidity after steam sterilization. The filters preserve their bacterial barrier, liquid barrier, and mechanical performance characteristics after sterilization.

The 9" x 9" disposable filter is designed for use with Genesis™ sterilization containers, sizes Full, Half, Mid, XL, and Retractor. They are equivalent to and can be used as a substitute for the Genesis DST-3 paper filter.

The filters of rigid reusable containers must cover the vented areas of sterilization ports to provide a microbial barrier and prevent recontamination of contents after sterilization.

Directions for Use:

- Place a filter over each vented area of the container
- Secure with the associated filter retention plate; lock in place
- Ensure that all perforated vents are covered
- Discard after each use

Caution:

- Filters are supplied non-sterile and for single use only
- Do not use more than one thickness of filter material

Storage Conditions

- Filters have no expiration date
- Keep clean and dry until use
- Store at room temperature away from moisture
- Keep covered within a plastic bag or in covered storage receptacle

Compatibility:

- Can be used wherever a 9"X9" polypropylene filter is required
- NOTE: Polypropylene filters are required for H₂O₂ sterilization cycles and may be used in Steam and EtO cycles



**SCF9S02
(500/PK)**