

SuperNova from Case Medical

Validated, controlled, effective 3 step cleaning process.

No Drums

No Wires

No Soap Closet



*Multi-Enzymatic Cleaner
CSNC25*

*Detergent
CSNB25*

*Lubricant
CSNL01*

Come out of the Closet with Case!

Save Money and Eliminate Storage Space

- ~ Precise Dosage Control
- ~ Efficacious 6 Log Reduction
- ~ pH Neutral Natural Products
- ~ Super Concentrated

Are you tired of moving 30 gallon drums? Try our SuperNova Cleaning System with Dose Control. Two of our 2.5 gallon totes equals 30 gallons of your current concentrated cleaner!

Introducing SuperNova[®], the world's most concentrated line of liquid detergents packaged in a 2.5 gallon tote with built-in handle. Highly effective for decontamination, and cleaning of medical devices. Multi-layered, pH neutral, enzymatic products and detergents can be used manually or in automatic washers to rapidly break down bio-burden and bio-film.

SuperNova



Validated at university laboratories for efficacy and biocompatibility. Our cleaning process demonstrated a ~6 Log reduction in bio-burden and >99.9% reduction of organic contaminants.



For manual cleaning, use SuperNova, one ounce equals 10 gallons of diluted product.



For mechanical cleaning, use SuperNova.25 formulated to meet the lowest default setting of most automatic washers.

SuperNova, a step-wise cleaning program utilizing the highest concentrated liquid enzymatic cleaners and detergents on the market followed by a highly concentrated instrument lubricant with rust inhibitor. Utilize our free dispensing system with proportioner when you purchase our line of SuperNova products.

***Please contact your Case Medical representative or customer service at
Case Medical, Inc.[®]***

65 Railroad Avenue • Ridgefield, NJ 07657

888-227-CASE • Fax: 201-313-9090

<http://www.casemed.com/>