

Case Medical

Universal Sterilization Containers

Case Medical has developed unique sterilization technology with the input of health-care professionals. The company's SteriTite sterilization container system has been validated for steam, gas, ozone and gas plasma/STERRAD Sterilization. Furthermore, the same SteriTite container can be easily adjusted for sealed flash sterilization by utilizing its optional FlashTite valve plate in lieu of disposable filter and retention plate. This product virtually eliminates the open pan for immediate sterilization of emergency items and full sets.

SteriTite containers are highly durable, heat treated, meet the MIL spec for corrosion resistance and are designed for aseptic presentation. Cost savings is built into the equation as our products are designed for more than one use and can reduce inventory through standardization. Plus, they are validated.

SteriTite is the world's first universal healthcare workers have a sealed container system that becomes a filter-less sealed flash container by merely exchanging the existing disposable filter and filter retention plate for the FlashTite

valve plate. The SteriTite container with FlashTite valve plate provides a significant improvement over other filter-less devices and virtually eliminates the hot, steamy open pan currently used for flashing near the OR. This increases sterility assurance as the instruments go into the autoclave sealed and are sealed within the container until point of use.

Another advantage of the SteriTite system is the unparalleled protection it affords when combined with the MediTray inserts. The MediTray products allow for versatility never before realized. No matter the set—whether endo, micro, orthopedic or Bookwalter—Case Medical has a validated solution. Loaner sets drop into our DIN sized case.

Customization is the added benefit in each of our offerings. The company's product line of sterilization trays includes both custom and standard components. Our Case Solutions can be customized for a specific use or application. Case Medical has patented and certified unique technology which places the company in the position of offering products that are safer, more effective and contribute to improved clinical outcomes.

For more information, visit <http://casemedical.com> or call 888.227.CASE.

