

SteriTite®

Case Carts and Shelving

SteriTite® case cart systems save time, money and effort. We offer open and closed case carts in single and double bay configurations. Designed for storage and transport, the universal size fits all sealed container systems. The stainless steel perforated shelves are fully adjustable and may be disassembled for cleaning. The shelf design avoids tears when removing wrapped sets from the cart or storage shelf.

Protect your instrument sets from contamination during transport and handling. All SteriTite carts and storage racks are manufactured from high quality stainless steel with adjustable perforated shelving, shock-absorbing wheels, cart washable casters. The flexible design maximizes your storage space and transportability. SteriTite case carts are designed to be cleaned manually and in automated cart washers. The durable construction and materials withstand high temperatures and chemical agents.

