

Introducing the All-in-One Mini-Pack



A full-scale reprocessing solution for smaller clinical settings

Provides the same free-rinsing, environmentally preferred, fast acting performance available in our standard sized products

Delivers a lower-cost strategy for achieving comprehensive results with validated cleaning solutions

Ideal for space constrained settings

U.S. EPA Safer Choice awarded formulation

Case Solutions are pH neutral instrument cleaners that meet AAMI/ANSI requirements for ideal cleaning agents. Inspired by customer requests for a pre-packaged assortment of cleaners, detergents, wipes, brushes and lubricants in sizes and quantities that are a better fit for clinical care practices, Case Medical is proud to offer you our latest kit.

Kit includes Safer Choice awarded cleaners:

- **SuperNova .25 Multi Enzymatic Cleaner** – 1 Gal
Ideal for long term soaking, manual and automated cleaning
- **SuperNova .25 Detergent** – 1 Gal
Non-enzymatic pH neutral cleaner for manual and automated cleaning
- **Penta Wipes** – 120 wipes
Ideal for surface cleaning and pre-cleaning of soiled devices
- **Alcohol Wipes** – 120 wipes
Ideal drying agent with some disinfection properties
- **PentaPrep Spray** – 24 oz.
Ready to use multi-enzymatic pre-cleaner validated for long term soaking
- **InstruCreme Lubricant** – 24 oz.
Aqueous emulsion compatible with devices, environmentally preferred
- **Micro/Eye Brush Kit**
In various sizes, color coded with blunted tip for lumens and narrow channels



Case Medical U.S. EPA Safer Choice Partner of the Year 2015, 2016 and 2017

© Copyright 2017 Case Medical, Inc.®