

SteriTite® Case Carts and Shelving

Designed of durable stainless steel to maximize internal space and hold containers for easy identification, SteriTite® Case Carts and Shelving are built for convenience and built to last. Shock-absorbing wheels and stainless steel bumpers help you get where you need to go with contents secure and intact.

Quality SteriTite® Case Carts

- **Quality Construction:** Durable stainless steel construction securely stores and transports surgical instruments and supplies
- **Configured for DIN-sized containers:** Depth accommodates containers for easy viewing of labels, indicators, and ID plaques
- **Available Disposition Monitor:** Permanently adhere to case cart to differentiate clean from soiled contents



CDS001 Disposition Monitor



SCCC01 Closed Single



SCCSD1 Closed Deep Wide

Versatile Open Carts and Shelving

- **Quality Construction:** Durable stainless steel construction transports and stores surgical instruments and supplies
- **Adjustable, Perforated Shelves:** Easy to clean and smooth, preventing holes or tears in blue wrap seen with wire shelves



SCC002 Open Double Bay Cart

SCC003H Open Triple Bay Shelving

Ordering Information - SteriTite® Case Carts and Shelving

Part ID	Description	Outer Dimensions w/ bumper - DxWxH	Outer Dimensions of shell - DxWxH	No. of Shelves	Shelves - DxW
Closed Case Carts					
SCCC01	Single	29"x22"x40"	26"x17"x33"	2	23"x12.6"
SCCC02	Double (2 bays)	29"x38"x40"	26"x33"x33"	4	23"x12.6"
SCCSD1	Short Deep	34"x33"x40"	30"x27"x33"	2	28.3"x23.3"
SCCTS1	Tall Standard	30"x33"x55"	26"x27"x48"	2	24"x23.3"
SCCTD1	Tall Deep	34"x33"x55"	30"x27"x48"	2	28.3"x23.3"
SCCC01CNB*	Under Counter	N/A	25"x21.5"x30.5"	1	17.6"x22.9"
Open Cart					
SCCO02	Double Bay Cart	27"x33"x44"	27"x33"x35"	4	26"x12.6"
Open Shelving					
SCCO03H	Triple Bay Shelving	27"x49"x65"	27"x48"x59"	9	26"x12.6"

*SCCC01CNB has no bumpers and is designed to fit under standard height counters



◀ **At Left:**
SteriTite
case carts
and shelving

At Right: SteriTite
Disposition
Monitors indicate
at-a-glance
whether contents
are **Clean** vs.
Contaminated
(order using Part ID
CDS001)



SteriTite is a registered trademark of Case Medical
Copyright © 2025 Case Medical Inc.
All rights reserved. Rev F 12/25